

DenTek Oral Care Products

Win a DenTek Bundle >>

Floss Picks + Case

Get a fresh clean mouth any time, anywhere.

With DenTek's Floss Picks + Case, flossing on the go has never been easier. Each sanitary case comes pre packed with 6 Comfort Clean Floss Picks and fits conveniently in pockets, purses or travel bags.

For the flosser on-the-go.

Don't get stranded with food stuck in your teeth!

- Have Germ Shield protection to prevent the growth of bacteria
 - Come in 5 vibrant colours
 - Are packed with 6 Comfort Clean Floss Picks per Case, which are great for tight teeth and sensitive gums
 - Are available in 2 case or 4 case packs



Instructions for Use

Flossing Instructions:

- Gently slide floss up and down between teeth allowing the floss to wrap around teeth.
 - Clean above and below the gum line.
 - Rinse as necessary.

Pick Instructions:

- Place pick between teeth and gently press against the gum.
- Use an in and out motion to remove food particles and stimulate gums.
- Rinse as necessary.

© 2014 DenTek Oral Care, Inc. All Rights Reserved.

v>