

DenTek Oral Care Products

Win a DenTek Bundle >>

Comfort Clean for Back Teeth Floss Picks

With Fluoride Coating

DenTek Comfort Clean for Back Teeth Floss Picks give you all the plaque fighting benefits of our regular Floss Picks, but are uniquely designed to easily reach the back teeth. With Soft Silky Floss, they are perfect for a comfortable clean between teeth.

DenTek Comfort Clean for Back Teeth Floss Picks come with all these benefits:

- Designed to easily reach back teeth
 - Silky Ribbon Floss
 - Fluoride coating
 - Mint flavour helps leave your mouth feeling fresh
 - Textured Pick cleans between teeth & helps stimulate gums
 - Tongue Cleaner helps to fight bad breath



Instructions for Use

Flossing Instructions:

- Gently slide floss up and down between teeth allowing the floss to wrap around teeth.
 - Clean above and below the gum line.
 - Rinse as necessary.

Pick Instructions:

- Place pick between teeth and gently press against the gum.
- Use an in and out motion to remove food particles and stimulate gums.
- Rinse as necessary.

Tongue Cleaning Instructions:

- Extend your tongue.
- Place the Tongue Scraper on the back of your tongue and pull forward, skimming the surface of your tongue.
- Rinse as necessary.

© 2014 DenTek Oral Care, Inc. All Rights Reserved.

v>