

# DenTek Oral Care Products

**Win a DenTek Bundle >>**

## DenTek Orabrush Tongue Cleaner

Helps fight bad breath

Up to 90% of bad breath comes from odour-generating residue on the top of your tongue. The DenTek Orabrush helps to remove this residue and keep your breath fresh.

Unlike toothbrush bristles which are designed for the smooth surfaces of your teeth, the DenTek Orabrush features ultra-soft, pointed bristles that reach deep into the uneven crevices of your tongue and loosen stinky bacteria. The scraper then collects and removes the bacteria and residue generated from brushing, helping to fight bad breath.

The long handle helps reduce gagging, and allows for an easier and more effective cleaning experience

Find out more about the Orabrush at: [www.orabrush.com](http://www.orabrush.com)



### Instructions for Use

- Moisten tongue with water or mouth wash.
- Gently brush tongue several times with the OraBrush using a back and forth motion reaching towards the back of the tongue.
- Remove the residue generated from brushing, by scraping your tongue with the scraper located at the top of

the brush.

- Rinse mouth to remove the loosened material leaving the tongue clean and fresh.
- Clean the OraBrush with water and store.

© 2014 DenTek Oral Care, Inc. All Rights Reserved.

v>