

DenTek Oral Care Products

Win a DenTek Bundle >>

Easy Brush ISO 1

Easily clean between tight spaces and dental work.

DenTek's Easy Brush ISO 1 makes cleaning between teeth and the tightest spaces easy. Flexible handle and wire make reaching trouble spots a breeze.

Reach trouble spots and dental work anywhere in your mouth. DenTek Easy Brush ISO1 reaches and cleans with:

- Durable tapered brush helps to clean tight spots; great for braces
 - Coated comfort wire bends for better access, but is easy on teeth and gums
 - Flexible Handle for comfort and control
 - Minty flavour that helps leave your mouth feeling fresh
 - Anti-Microbial cap helps keep brush clean between uses
 - Cap fits on grip to extend the handle's reach



Instructions for Use

- Insert between teeth, interdental spaces, or dental appliances.
- Use a soft back-and-forth motion to remove plaque and food particles.
- Rinse after use.
- Replace when bristles become worn.

© 2014 DenTek Oral Care, Inc. All Rights Reserved.

v>