

# DenTek Oral Care Products

**Win a DenTek Bundle** >>

## Easy Brush ISO 2

Easily clean between spaces, braces, and dental work.

DenTek Easy Brush ISO 2 cleans between teeth and hard to reach trouble spots. Great for cleaning around braces, bridges, and other dental work.

Reach trouble spots and dental work anywhere in your mouth. DenTek Easy Brush ISO 2 reaches and cleans with:

- Durable tapered brush helps to clean tight spots; great for braces
  - Coated comfort wire bends for better access but is easy on teeth and gums
  - Flexible Handle for comfort and control
  - Minty flavour that helps leave your mouth feeling fresh
  - Anti-Microbial cap helps keep brush clean between uses
  - Cap fits on grip to extend the handle's reach



## Instructions for Use

- Insert between teeth, interdental spaces, or dental appliances.
- Use a soft back-and-forth motion to remove plaque and food particles.
- Rinse after use.
- Replace when bristles become worn.

© 2014 DenTek Oral Care, Inc. All Rights Reserved.

v>