

DenTek Oral Care Products

Win a DenTek Bundle >>

Fun Flossers

Make flossing fun!

DenTek Fun Flossers make flossing fun for children while teaching healthy oral care habits for life.

All DenTek Floss Picks help reduce tooth decay and remove food and plaque.

DenTek Fun Flossers encourage flossing:

- Easy to grasp handle for both kids and parents
- Multi-filament textured floss fits kid's teeth
- Small flosser head fits children's mouths
- Wild Fruit flavour that kids love



Instructions for Use

- Hold flosser firmly between thumb and index finger.
- Using a gentle back-and-forth motion, slide floss between teeth.
- Move up & down against tooth to clean above & below the gum line.

- Repeat steps for each tooth.
- Discard after use.

Warnings:

Gums may bleed when you first begin to floss. If bleeding persists or gums remain red and swollen, see your dentist. Adult supervision of children is necessary to ensure that they floss correctly.

© 2014 DenTek Oral Care, Inc. All Rights Reserved.

v>