

DenTek Oral Care Products

Win a DenTek Bundle >>

Sensitive Extra Gentle Floss Picks

With Fluoride Coating

DenTek Sensitive Extra Gentle Floss Picks give you all the plaque fighting benefits of our regular Triple Clean Floss Picks, but are designed with a Soft Silky Ribbon Floss, perfect for sensitive teeth and gums.

- Soft Silky Ribbon Floss with Fluoride
- Easy on Sensitive Teeth & Gums
- Textured Pick cleans between teeth & stimulate gums
- Tongue Cleaner helps fights bad breath with a minty flavour



Instructions for Use

Flossing Instructions:

- Gently slide floss up and down between teeth allowing the floss to wrap around teeth.
 - Clean above and below the gum line.
 - Rinse as necessary.

Pick Instructions:

- Place pick between teeth and gently press against the gum.
 - Use an in and out motion to remove food particles and stimulate gums.
 - Rinse as necessary.

Tongue Cleaning Instructions:

- Extend your tongue.
- Place the Tongue Scraper on the back of your tongue and pull forward, skimming the surface of your tongue.
- Rinse as necessary.

© 2014 DenTek Oral Care, Inc. All Rights Reserved.

v>