

DenTek Oral Care Products

Win a DenTek Bundle >>

Tongue Cleaner

Don't surrender to bad breath.

DenTek's Tongue Cleaner helps to scrape away bad breath odours at the source.

Fresh mint flavour.

Fresh breath made easy. DenTek Comfort Clean Tongue Cleaner helps to scrape away bad breath odours at the source. And its low profile and narrow design comfortably reach the back of the tongue, so you don't have to gag to fight bad breath.

- Effective way to remove bad breath causing bacteria from hard to reach places
- 3 cleaning edges ergonomically designed for cleaning the curves of the tongue



Instructions for Use

- Extend your tongue.
- Place the Tongue Cleaner on the back of your tongue and pull forward, skimming the surface of your tongue.
- Rinse the Tongue Cleaner with warm water and repeat as necessary.
- If bad breath persists, consult your dentist.

For Best Results:

- Like your toothbrush, you should replace your tongue cleaner every 3 months.

© 2014 DenTek Oral Care, Inc. All Rights Reserved.

v>