

DenTek Oral Care Products

Win a DenTek Bundle >>

Triple Clean Floss Picks

Three ways to get your teeth & mouth cleaner; floss, pick and tongue cleaner.

DenTek Triple Clean Floss Picks gives you three ways to get your mouth cleaner, floss, pick and tongue cleaner. Shred resistant textured floss is ready to conquer the tightest of teeth.

Want a floss pick that scrubs between tight teeth, helps stimulate gums, and helps remove the germs that cause bad breath better than any other floss pick?

DenTek Triple Clean is the floss pick for you. Get three times the clean:

- Shred Resistant Textured Floss
 - Great for tight teeth
 - Textured Pick cleans between teeth and helps stimulate gums
 - Tongue Cleaner helps fight bad breath with a minty flavour



Instructions for Use

Flossing Instructions:

- Gently slide floss up and down between teeth allowing the floss to wrap around teeth.
 - Clean above and below the gum line.
 - Rinse as necessary.

Pick Instructions:

- Place pick between teeth and gently press against the gum.
- Use an in and out motion to remove food particles and stimulate gums.
- Rinse as necessary.

Tongue Cleaning Instructions:

- Extend your tongue.
- Place the Tongue Scraper on the back of your tongue and pull forward, skimming the surface of your tongue.
- Rinse as necessary.

© 2014 DenTek Oral Care, Inc. All Rights Reserved.

v>