

DenTek Oral Care Products

Win a DenTek Bundle >>

Whitening Floss Picks

With Fluoride Coating

DenTek Whitening Floss Picks give you all the plaque fighting benefits of our regular Floss Picks, but are designed with Whitening Micro Crystals and Soft Silky Floss, perfect for whitening between teeth.

DenTek Whitening Floss Picks come with all these benefits:

- Whitening Micro Crystals
- Soft Silky Ribbon Floss
- Fluoride coating
- Mint flavour helps leave your mouth feeling fresh
- Textured Pick cleans between teeth & helps stimulate gums
- Tongue Cleaner helps to fight bad breath



Instructions for Use

Flossing Instructions:

- Gently slide floss up and down between teeth allowing the floss to wrap around teeth.
- Clean above and below the gum line.
- Rinse as necessary.

Pick Instructions:

- Place pick between teeth and gently press against the gum.
- Use an in and out motion to remove food particles and stimulate gums.
- Rinse as necessary.

Tongue Cleaning Instructions:

- Extend your tongue.
- Place the Tongue Scraper on the back of your tongue and pull forward, skimming the surface of your tongue.
- Rinse as necessary.

© 2014 DenTek Oral Care, Inc. All Rights Reserved.

v>